










Ball Attachment

The ball attachment superstructure is intended to secure a removable prosthesis. The attachment is used in conjunction with a stainless steel cap and an intermediate silicone insert.





Material: Titanium (Ti6Al4V ELI)



Ball Attachment	L (mm)
 NM-T1200	0.5
 NM-T1201	1
 NM-T1202	2
 NM-T1203	3

Ball Attachment	L (mm)
 NM-T1204	4
 NM-T1205	5
 NM-T1206	6

Caps for 2.5mm Ball Attachment

Metal	Extr Soft Silicon	Soft Silicon	Standard Silicon
 NM-T3001	 NM-T3002	 NM-T3003	 NM-T3004